



Release Notes for the BR.io Cloud Platform

Version 1.0

October 2020

Contents

Introduction.....	1
Supported Operating System.....	1
Summary of Main Features	2
Known Issues	4
Documentation	4
Legal Notices.....	5
Contacting Technical Support.....	5

Introduction

Using the BR.io Cloud Platform with CFX Opus instruments you can set up, run, monitor, and analyze your real-time PCR experiments on connected CFX Opus instruments. This document summarizes the main features, and also provides a short list of known issues.

Supported Operating Systems and Browsers

BR.io is supported on

- Windows 7 OS or later
- Mac OS 10.12 or later
- Chrome browser v63 or later
- Safari browser v11 or later

Summary of Main Features

Individual User Accounts

- You can create an individual BR.io user account, remain logged in for up to 30 days, change a password, and reset a forgotten password.

Online Assistance and Accessibility

- BR.io features a Help Center that provides searchable software documentation and context-sensitive Help that is coordinated with primary BR.io pages.
- The BR.io interface is designed for ease-of-use by all users.
- Users can contact Bio-Rad Customer Support directly from BR.io.

Storing and Managing Data

- BR.io displays the full list of a user's files on the Files page, and the most recently used files on the Home page.
- Users can upload .zpcr and .pcrd file types to the Files page, and rename and delete files.

Integration with CFX Opus Instruments

- Using the BR.io Cloud Platform you can
 - Connect CFX Opus instruments to your BR.io account
 - Create a run in BR.io, and then execute the run on the CFX Opus instrument
 - Create a run on the CFX Opus instrument using a locally stored protocol, and save it to your storage area in BR.io
- For runs created or saved in BR.io, CFX Opus run data is automatically uploaded to BR.io from connected instruments.

Managing and Monitoring CFX Opus Instruments

- In BR.io, you can view
 - A list of all CFX Opus instruments connected to your BR.io account
 - The model, name, and serial number of a connected CFX Opus instrument
 - The current status, time remaining, and current step/cycle for an experiment running on a connected CFX Opus instrument

Creating and Managing CFX Protocols

- BR.io provides a dedicated repository for your protocols.
- You can create, edit, and delete protocols in BR.io.
- CFX protocols in BR.io appear in both a graphical step format and a list of protocol steps.
- In BR.io, you can create temperature, gradient, melt curve, and goto steps in your protocols.
- BR.io validates your protocol parameters and displays an error message if applicable.

Creating and Editing a CFX Opus Run in BR.io

- BR.io provides a
 - Streamlined workflow to set up an experiment to run on the CFX Opus 96 or CFX Opus 384
 - Plate editor that is optimized for usability, where you can see well content, select and edit multiple wells at once, and see plate validation and error messages, if applicable

Note: BR.io only supports the following scan modes:

SYBR/FAM only
All channels

- You can edit the plate and run details while the run is in progress, or after it has completed on the instrument.
- You can set up a CFX Opus run using a default (quick) plate. Default plates consist of pre-defined fluorophores and unknowns in every well, so you can skip or defer plate setup.

BR.io Analysis Features

- BR.io uses the same algorithms as CFX Maestro.
- BR.io enables the user to analyze real-time PCR (amplification) data, including the ability to
 - Toggle baseline subtraction on or off
 - Set automatic or user defined single thresholds
 - View and download the amplification chart and apply linear or logarithmic scaling of the RFU data
 - View amplification data in a sortable table format
 - Export amplification data table and RFU data per cycle in .csv format
 - Analyze data by fluorophore or target
 - Show or hide data from individual wells in the plate
 - Filter out targets or fluorophores from analysis
- Analyze melt curve data from CFX runs (Opus and legacy), including
 - A melt peak (negative derivative) chart that you can download and print
 - A sortable, exportable data table, which displays melting temperature, peak height, begin temperature, and end temperature

Interoperability with CFX Maestro or CFX Manager

- You can
 - Upload and import .zpcr files and .pcrd files created in CFX Manager v1.0 and later, or CFX Maestro v1.0 and later
 - Upload and import .pcrd files from CFX Maestro for Mac
 - Download .zpcr file from a CFX Opus run
- Note:** The plate layout and run details are omitted.
- Download the original .zpcr or .pcrd file from a manually uploaded (imported) CFX run

Known Issues


- Users who mistype their email address during signup can experience issues with signing in or resetting their password.
- You must close the “Run successfully uploaded to your BR.io account” dialog box soon after the run is completed, or BR.io incorrectly displays the CFX Opus status as Offline.
- If you navigate from the CFX run workflow while uploading a file, BR.io does not warn you about unsaved changes to your CFX run.
- After you exit the Target Name field, BR.io assigns a color to a new target on the Plate Setup for a CFX run.
- Plate headers (row and column) do not remain in view on Plate Setup while scrolling.
- BR.io does not currently support
 - Analysis of .pcrd or .zpcr files that contain legacy or user-calibrated fluorophores
Note: You can upload the files, but working with them in the Analysis module can produce errors.
 - Viewing or analyzing completed CFX runs with melt curve data, but no amplification step
 - Viewing the gradient (temperature per row) when adding a Gradient step to a CFX protocol
 - Viewing details and metadata for a completed CFX run, including the date and time of the run, and the name, model and serial number of the instrument
 - Downloading CFX runs created or modified in BR.io as CFX Maestro-compatible PCRD files
 - Baseline adjustment and fluorescence drift correction analysis settings
 - Application-based analysis, such as standard curve/absolute quantification, gene expression/relative quantification, and allelic discrimination

Documentation

Click the following URL, and then click the  icon to access the online Help Center.

<https://br.io>

You can also click the down arrow above EXPLORE BR.io to access videos that describe BR.io processes.



Legal Notices

No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording, or any information storage or retrieval system, without permission in writing from Bio-Rad Laboratories, Inc.

Bio-Rad reserves the right to modify its products and services at any time. This guide is subject to change without notice. Although prepared to ensure accuracy, Bio-Rad assumes no liability for errors or omissions, or for any damage resulting from the application or use of this information.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

SYBR is a trademark of Thermo Fisher Scientific, Inc. Bio-Rad Laboratories, Inc. is licensed by Thermo Fisher Scientific, Inc. to sell reagents containing SYBR Green I for use in real-time PCR, for research purposes only.

All trademarks used herein are the property of their respective owner.

Copyright © 2020 Bio-Rad Laboratories, Inc.

Contacting Technical Support

The Bio-Rad Technical Support department in the U.S. is open Monday through Friday, 5:00 AM to 5:00 PM, Pacific time.

Phone: 1-800-424-6723, option 2

Email: Support@bio-rad.com (U.S./Canada Only)

For technical assistance outside the U.S. and Canada, contact your local technical support office or click the Contact Us link at www.bio-rad.com.